# REACTION Software Ver.9 PC Requirements



## Is a New Computer/PC Necessary?

Today's PG1000 REACTION software requires a high-performance PC to run properly. Older versions of our software captured images and applied single calculations to the images whereas today's software is taking hundreds of calculations per second on live images. This requires a significant amount of computing power.

Computer details can be found under System Information or by right clicking on My Computer or This PC and selecting properties.

- Operating System
  - o 32 bit must replace PC
  - o 64 bit OK
- Processor
  - o i3 must replace PC
  - i5 upgrade suggested, performance issues (i.e. lagging) may occur. Crashing may occur when using edge detection features.
  - **i7 OK**
- RAM
  - o 2GB must replace PC
  - 4GB upgrade suggested, performance issues (i.e. lagging) may occur. Crashing may occur when using edge detection features.
  - 8+GB OK
- Windows
  - Windows 8/Vista/XP must upgrade to Windows 10
  - Windows 7 must upgrade to Windows 10
  - Windows 10 OK
  - Windows 11 NOT TESTED
- Graphics Card
  - o Required for 4K systems/model UI358x cameras

We can not upgrade software during a service call if the PC does not meet these minimum requirements.

### **Upgrade Information**

Part No.	Description
00087-5P	High-performance PC with REACTION Software preloaded

#### Contact PG Inspection Technologies for a quotation or for additional help.

Please have your PG1000 model and serial number, camera model and serial number, software version number and computer information available.

www.pg1000.com

pg@pg1000.com

(262) 946-5420

#### **Examples:**

System

Manufacturer: Dell

Model: Optiplex 390

Rating: 5.0 Windows Experience Index

Processor: Intel(R) Core(TM) i3-2120 CPU @ 3.30GHz 3.30 GHz

Installed memory (RAM): 3.00 GB

System type: 32-bit Operating System

This computer will need to be replaced. The current REACTION software will not run on 32-bit operating systems.

Device Manager Windows edition Remote settings Windows 7 Professional System protection Copyright © 2009 Microsoft Corporation. All rights reserved. Advanced system settings Service Pack 1 Get more features with a new edition of Windows 7 System Manufacturer: Dell Model: Optiplex 3010 Rating: 5.3 Windows Experience Index Processor: Intel(R) Core(TM) i5-3450 CPU @ 3,10GHz Installed memory (RAM): 4.00 GB System type: 64-bit Operating System

This computer will require a Windows 10 update but even after the update there may be performance issues (image lag or crashing) due to the older i5 processor and only having 4GB of RAM. Ideally the user would upgrade to an i7 processor, 8-16GB of RAM and Windows 10Pro.

# View basic information about your computer

Windows edition

Windows 10 Pro

© 2019 Microsoft Corporation. All rights reserved.



System

Processor: Intel(R) Core(TM) i5-8500T CPU @ 2.10GHz 2.11 GHz

Installed memory (RAM): 8.00 GB (7.78 GB usable)

System type: 64-bit Operating System, x64-based processor

This computer should work. It has a newer i5 processor and 8GB of RAM as well as Windows 10.